

(For the use in Medical profession)

**Anmol Pulsar V - 5
(Kumar's Magnet)**

Customer care: 9247110759


V5 for Women's Health

Most important instructions: do not use during:

- 1. Once conceived / or a period is missed*
- 2. During Pregnancy*
- 3. Two days before, during and after Menstruation*
- 4. Any other situation which you find doubtful*

This is a Magnet therapy device useful in Physiotherapy. To be used under the guidance of an experienced Magnet therapist or Physiotherapist or a Doctor.

This is not intended for primary medical care but only as a supportive one. Consult your doctor before using. Most of the models are specifically designed for individual usage, hence be careful while you get one for yourself.



This Anmol Pulsar V-5 is designed to cater the needs of Women's health. This can be used along with other / conventional therapy / treatment. The synergistic effect is observed to last more and comes in a rapid pace hence reducing the intake of medicines.

The invisible and non-invasive magnet therapy based on PEMF principles gives better results.

Anmol Pulsar V-5 is the latest generation of Magnet therapy devices with PEMF (Pulsating Electro Magnetic Field) principles which has various power sources to give better relief at most times and places.

Anmol Pulsar V-5 is easy to use; no much prior training is required other than some mandatory cautions to be followed. Can be used for most of the people from ages 12 (twelve) and above.

Anmol Pulsar V-5 is dual powered with solar and grid (household) power. However it's much recommended to use on solar as the overall efficacy is found to be much higher than when compared to the same tenure of exposure taken at night or with grid power.

Anmol Pulsar V-5 works on Pulsating Electro Magnetic Field principles and brings the effect on the area it engulfs of the subject to treatment. The magnetic field only engulfs the area but the current is not induced forcibly into the body, hence it's much safer than other conventional TENS or similar type of devices. Though there are limitations for every device or concept, this Anmol Pulsar V-5 surely is the less risk free option when handled by inexperienced or less knowledgeable persons.

All one need to do is to place the device's probe / plate near the affected area for a few minutes daily. The solar panel has to be placed in a well sun exposed area and the cable drawn to the place of treatment.

Anmol Pulsar V-5 helps to:

1. Absorb the externally applied medication on the skin for better action
2. Relieve acute pains
3. Improves health and immunity
4. Safe, easy, simple to handle or apply (self operative)

A must for every home, clinic, hospital and a wise asset.

Precautions:

1. Should be used under the guidance or in accordance with your family doctor or with a person who knows how to use this Anmol Pulsar V-1 safely and with cautions.
2. Should not be used during pregnancy
3. Should not be used for children below 12 years
4. Should not be used with persons having electronic implants
5. Should not be used with persons having braces and surgical implants
6. Should not be used during menstrual period, TB, viral diseases, internal hemorrhage, fevers, heart ailments etc.
7. Should not be placed directly on / near head as this is to be applied on the navel / pubis area for women's health.
8. Should not be used as a replacement for conventional medical treatment or as a reason to postpone seeing your doctor.

9. Should not be used when there is bleeding or on an undressed injury.
10. Consumption of magnetic water too should be limited to maximum of three times / day
11. Take all caution to place the head towards south or west sides ONLY while taking position.
12. If any sweating or irritation or any uncomfortable situation arises discontinue the usage and consult a doctor
13. A maximum of 20minutes exposure for a maximum of three times a day is necessary but in most cases just a 20 minute exposure once in a day is enough followed with same time of undisturbed movement of the body for the maximum effect to take place.
14. In most situations the blood circulation improves and people may get sleep during exposure / treatment, it is recommended not to disturb them at those times.
15. Taking of hot liquids is recommended before and after exposure / treatment. Taking cold liquids / beverages SHOULD BE AVOIDED before / after treatment / exposure.
16. There should not be any electronic items like mobiles, tablets, computers etc within a vicinity of about 3 feet from the device. Since this Anmol Pulsar V-1 radiates low frequency emissions, it may harm other devices.
17. Keep the equipment out of reach of children
18. This unit contains a powerful magnet; hence see that any electronic or electromechanical devices are away from the device. They might get damaged permanently
19. If you are not sure where or how to apply contact customer care or contact your family doctor or contact anyone well versed with magnet therapy.

Effective on:

Issues related to lower part of body from navel to toes, Back ache, knee joint pains, burns, wounds, ulcers, women's issues and other circulatory disorders.

Specialties:

Improves blood circulation, health and immunity, helps to normalize disorders of nervous and vascular system. Improves natural perspiration of skin.

Merits:

Safe replacement for hot packs. No warming of tissues. No needles used. No need for primary electrical impulses. Non-hazardous and no unpleasantness. No need for specially trained persons to handle.

Model: Anmol Pulsar V5 Specifications:

1. Solar panel of 12volts 4watts – 1no (model without solar panel is also available)
2. AC adaptor of 12volts – 1no
3. Energy device convertor – 1no
4. Probe – 1no
5. High power magnets in separate pouches – 2nos
6. High power magnets for foot application – 2nos


Fixing and using the unit:

If you are using the unit on solar mode, see that the solar panel is fixed in a well lit and sun exposed area and the solar module / panel is aptly fixed so that air / wind issues will not disrupt the panel position. See that the solar panel is fixed with a little inclination so that in the event of rain, fog etc they will slide down easily. The best direction to fix the panel is

south facing and inclined. Then run the wire / cable attached to the panel to the place where the subject is in place. The solar panel cable has a pin, attach this to the converter unit. Once you connect the cable switch on the device. You should now see a light Green will glow, if the probe is properly functioning you can feel a little vibration, which indicates that the unit is ready and functioning.

If you want to run the unit with household power with the adaptor, connect the adaptor to the wall socket and switch on the AC power. Connect this cable to the energy converter and switch on the device. You should now see a light Green glowing, if the probe is properly functioning you can feel a little vibration, which indicates that the unit is ready and functioning.

Let the subject / person lie down on supine position & head towards south or west direction.



Placement of Probe and Magnets:

First set the power connections and arrange the unit after carefully integrating all the components. Keep the palms and feet separately, don't keep them crossed. Keep the head towards south or west sides only, so that you will see north or east. Treatment time is 20minutes for any session or as recommended. **2 glasses of water should be taken before treatment.** Avoid taking the treatment immediately after food. A gap of 30 min is recommended before and after food. Do not place the legs directly on the floor immediately after treatment, wait for 10 minutes and then use slippers / footwear for the next half an hour. Cold food and drinks should not be consumed for the next 30min or so.

Place the probe between the navel and pubis. Keep the two square magnets below the knees. Tie the foot magnets to the feet, one inside the ankle and other outside the ankle.

Duration of exposure / treatment:

A maximum of 20 minute exposure is enough per sitting with a maximum of three sittings per day in equal intervals. In case of joint pains use of external medicated oils or any sort of rubefacient are an added advantage. See that the precautions are met without fail. Though the unit is water resistant (not water proof) see that no part comes in contact with water or any other sort of liquid. All Anmol Pulsar Versions does not come with any controls hence there is nothing to manipulate and is designed for care as per the versions specifications.

For uterine related issues:

1st week: Place the probe on navel / pubis and square magnets below the back, near to kidneys and small magnets below the knee, once or twice a day as required.

2nd week: Place the probe on navel / pubis and square magnets below the knee and small magnets to the bottom of foot, once or twice a day as required.

3rd week: Place the probe on navel / pubis and square magnet below the ankle, tie the small magnets one on inside the ankle and the other outside the ankle, once or twice a day as required.

If the condition still persists, repeat the cycle. Also continuing with Ayurvedic medicines is highly recommended as both go hand in hand for best synergistic effect.

For digestive, kidney, hip issues:

Morning session: Place the V5 Probe on the navel (while in supine position) and place both the magnets on either side of the back / behind the kidney areas.

For knee joint pains and sciatica:

Morning session: V5P2 – Vertically place the V5 Probe (while in supine position) in between the knees and both the magnets on either side (outside of) knees, vertically. Applying Anmol oil is complementary.

Or

Morning session: V5P2 – Place the V5 Probe on the navel (while in supine position) and both the magnets below the knees. Applying Anmol oil is complementary.

Evening session: V5P1 – Place the V5 Probe on the navel (while in supine position) and both the magnets below the back / near kidneys. Applying Anmol oil is complementary.

For ankle joint pains and other issues in that area:

Morning session: V5P3 – Vertically place the V5 Probe (while in supine position) between the ankles and both the magnets on either side (outside of) ankles. Applying Anmol oil is complementary.

Evening session: V5P1 – Place the V1 Probe on the navel (while in supine position) and both the magnets below the back / near kidneys. Applying Anmol oil is complementary.

When to use:

1. It is always better to take treatment / exposure before one hour of intake of food.
2. Try to maintain the same time everyday for better efficacy
3. Once a while treatment is not effective
4. In case of only for rejuvenation use for 20minutes a day either sitting on it (on the probe) with necessary cushion.
5. For elderly people the probe can be placed below the foot for 20minutes a day, daily
6. Should not be applied near to head and heart, this is not designed to handle those issues.
7. After the sitting of 20minutes, its advised to take rest for another 20minutes without any movement
8. The best position to use is to place the probe below the navel point in the supine position

Observation during treatment:

1. There might be some drawing sensation in limbs after 5-10minutes of the commencement of treatment. This happens only in initial stages that too only for some. If it's unpleasant withdraw the sitting and consult your doctor.
2. As the blood circulation starts improving, the treated may experience slight palpitation of heart. If this happens for a longer time discontinue the treatment and consult your doctor. The feeling of body becoming light and tendency to sleep may occur. For many the symptom of sleep continues.
3. After the sitting, the treated will feel rejuvenated and experiences activeness in routine life.

Magnetized water:

After getting the desired results, it is better to drink magnetized water at least once a day for a minimum of 15 to 20 days or every day ever after. Keep one or two ounces of fresh drinking water filled in a glass container on the probe for 10 to 20 min and the same can be consumed. The probe should be switched on. This magnetized water is very much helpful for proper maintenance of health. Do not store this water for a long time. It is to be consumed in 30minutes.

Magnetized water can also be used for household plants. This is particularly useful for creepers and crawlers and flowering plants.

Troubleshooting:

1. If the probe is not vibrating check the pin of the power input of the power converter. If the connections are okay then the power converter indicator too will light up or else the indicator will not light up.
2. If the indicator on the power converter light up and only the probe is not vibrating then refer the unit to the supplier
3. If using on solar panel mode see that the best time to operate is between 10am to 4pm
4. If you find any of the units (probe or the power converter) are hot which is an anomaly, refer to servicing
5. Do not use the adaptor given to any other device or for any other purpose.
6. Do not use the solar panel for other necessities.

Notice: This device is patented and all the instructions are given as per the findings from individual research and general data. Any infringement and usage of information provided and if the designs are copied will lead to legal complications. Use of this device is at your own advantage and discretion. To be used only upon advice of your family doctor or a registered medical practitioner or magnet-therapist.

Sold and serviced by:

Sri Sresti Marketing

16-3-32/B, Kowtha Street,
Purnanandampet,
VIJAYAWADA – 520003. Andhra Pradesh.
Phone: 927110759

Email: srestimarketing@gmail.com website: www.srisms.com